

BioParaFree

 Catalog number
 Unit size

 08-1750
 100 ml

 08-1750-X4
 4 x 100 ml

Packaging Primary container: White bottle in PEHD with atomizer.

Secondary container: carton box.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Expected aim Product for the preparation of cyto-histological samples for optical microscopy.

Specifications Density: 848 Kg/m3 (20°C)

Viscosity: 15 cSt (40°C) Flash point 190°C

Insoluble in water and soluble organic solvents.

Application Microtomys and workbenches cleaning. It used for the removal of the residual paraffin.

Method 1) Atomize on the cleaning area, after 10 - 20 seconds, remove the residual paraffin.

Components CAS CE Index

Mixture of mineral oils

Functional properties

Mixture of mineral oils of high degree of purity according to FDA quality requirements.

The product is odorless and colorless.

Warning and precaution

The product must be used exclusively by specialized technical operators.

Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental

release.

Do not use if the primary container is damaged.

Storage Store in a cool dry area, at room temperature. Keep away from sources of ignition.

Keep from light. Keep the containers tightly closed.

Stability After the first opening, the product is usable until the expiry date, if correctly stored. Validity: 2 years.

Disposal Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

Date of issue: July 2013

Producer: Bio-Optica Milano s.p.a. Page 1 di 1